

1Cell.AI Oncology Solutions

Comprehensive Genomic Testing Portfolio

This booklet presents 1Cell.AI's comprehensive portfolio of AI-powered genomic testing solutions designed specifically for oncologists managing complex, advanced, and refractory cancer patients. Each test is built on cutting-edge NGS technology and validated clinical evidence to provide actionable insights that drive precision-based treatment decisions.

What we offer

Comprehensive Tumor Profiling	Liquid Biopsy Solutions	Disease-Specific Advanced Panels	Monitoring & Recurrence Assessment
Deep molecular insights from tissue samples to guide therapy.	Non-invasive testing for real-time tumor assessment.	Focused solutions for specific cancers and clinical questions	Track disease status over time.

Test Portfolio Overview

Test	Focus Area	Key Features
OncolIdx	Comprehensive genomic profiling	1000+ gene panel, actionable mutation detection, comprehensive analysis for solid tumors
OncolIdx Prime+	Tissue-Liquid-Normal matching	Tissue-liquid-normal matching, expanded biomarkers, supports rare cancers & CUP
OncolIdx TBx	Integrated genomics + transcriptomics	DNA + RNA analysis, fusion detection, integrated insights
OncoTarget	Targeted disease-specific tests	Customizable panels, rapid and cost-effective, disease focused approach
OncoHRD	Homologous Recombination Deficiency	Detects HRD, analyzes 50+ HRR genes, supports PARP inhibitor therapy
OncoRisk	Hereditary cancer risk assessment	Germline testing, 74-gene panel, familial cancer screening
OncoMonitor	Longitudinal cancer monitoring	Liquid biopsy, MRD detection, tracks, treatment response tracking, relapse prediction

Key Benefits of Our Test Portfolio

Precision Medicine: AI-enhanced analysis and comprehensive gene panels provide clinically actionable insights tailored to each patient's tumor profile.

Validated Science: All tests are built on rigorous clinical and technical validation with published performance characteristics.


Comprehensive Coverage: From initial diagnosis through treatment monitoring, our portfolio supports the entire cancer patient journey

Clinical Integration: Reports designed for easy interpretation by oncologists, with evidence-based therapeutic recommendations.

Product Portfolio supporting lifecycle of patient

Supporting You Across the Cancer Care Journey

Our portfolio is designed to ensure that at every critical decision point, you have access to the right test, at the right time, for the right patient.



Beyond Testing:
Clinical Decision Support

Dr G

A dedicated clinician app designed to simplify report interpretation.

MTB

Access to a multidisciplinary panel of experts to interpret complex cases

OncoIndx[®] offers comprehensive genomic profiling across solid tumors. Its unique design identifies key actionable insights for personalising cancer treatment



OncoIndx[®] in practice

- Challenges of managing complex cancers are addressed through evidence-based treatment choice.
- NGS report enhanced with artificial intelligence and updated clinical studies help identify the most beneficial therapeutic options.
- An easy digital interface to interpret and understand complex cancer genomics.

Where is OncoIndx[®] needed?

- Cancers progressing on standard-of-care treatment.
- Cancers where personalised treatment could enhance survival.
- Cancers showing resistance to conventional treatment.
- Cancers where combination of target drugs and/or immunotherapy could benefit.

Clinical and Technical Validation of OncoIndx[®] Assay A Comprehensive Genome Profiling Assay for Pan-Cancer investigations

Alteration	* PPV	* NPV	Accuracy	Specificity	Sensitivity
SNVs	100	100	100	100	100
Small INDELs	100	94.03	97.40	100	95.60
CNA	100	100	100	100	100
Fusions	100	96.43	98.48	100	97.44


* NPV: Negative predictive value. PPV: Positive predictive value.

Scan QR to read
the full article



Designed for oncologists managing complex, advanced and refractory cancer patients

- Proprietary panel of 1,000+ genes
- Longitudinal monitoring markers
- PD-L1 through IHC (TBx)
- CTC enumeration
- Traceability to ongoing clinical trials worldwide



Key Actionable Genes Include:



Lung
EGFR, ALK, HRAS, NRAS, KRAS, BRAF, TERT, ROS1, MET, RET, NTRK1/2/3, ERBB2, PIK3CA



Melanoma
BRAF, KIT, ROS1, NTRK1/2/3



Colorectal
BRAF, KRAS, NRAS, HRAS, ERBB2, MLH1, MSH2, MSH6, PMS2, NTRK1/2/3



Gastrointestinal
ALK, HRAS, HER2, NRAS, KRAS, BRAF, TERT, ROS1, MET, RET, NTRK1/2/3



Prostate
BRCA1, BRCA2, MLH1, MSH2, MSH6, PMS2, AR, NTRK1/2/3



Breast
BRCA1, BRCA2, HER2, PIK3CA, MSH6, PMS2, NTRK1/2/3, MLH1, BRAF



Bladder
FGFR2/3, NTRK1/2/3



Ovary
BRCA1, BRCA2, HER2, PMS2, MLH1, MSH2, MSH6, NTRK1/2/3



Thyroid
ALK, HRAS, NRAS, KRAS, BRAF, TERT, RET, NTRK1/2/3

Features of Oncolndx report

- In accordance with ACMG/AMP/ASCO/CAP guidelines
- Clear, unambiguous and actionable genomic findings.
- Based on highly reliable and internally curated databases.

Technical Specification

Specimen	Liquid	Tissue
Panel Size	3 MB	3 MB
Sample Type	cfDNA (TagFlexLB/Norgen/PaxGene Tube) + 1 tube of 8.5 ml Peripheral Whole Blood (~8-10 days of washout period required if patient under therapy e.g., Chemotherapy, Radiation, Immunotherapy, etc.) Peripheral blood using butterfly needles/ 22 gauge syringe, CSF, Plasma, Pleural Fluid	1 FFPE tissue block with 20-30 % tumor content (not more than 3-4 months old)* FFPE tissue block
Sample volume/quantity	8.5 ml (2 tubes)	1 FFPE block with 20-30% tumor content
cfDNA/DNA input amount	> 1ng	> 1ng
Sequencing platform	Illumina	Illumina
Read Length	2 x 150bp	2 x 150bp
Mean Sequencing Depth	10000x	2000x
Sequencing Chemistry	Hybrid capture	Hybrid capture
Indications	Solid Tumor	Solid Tumor
Analysis and Reporting	iCare™- AI based analysis platform with active reporting options	

Oncolndx® is a well-validated NGS test to extract therapy information from tumor DNA. It is possible on both tissue as well as blood sample.

1000+ proprietary genes panel

This NGS test includes 1000+ genes of targetable/ therapeutic significance.

>99% on target coverage

The test has a mean coverage of 98-100% across all tested genes.

2000x / 10000x mean depth

The test is performed at a mean depth of 2000x for tissue samples and 10000x for liquid.

OncoIndx[®] Prime⁺ is a Comprehensive multi-technology, multi-specimen test which performs molecular profiling of tumor and normal samples simultaneously from the same patient



Why Tissue-Liquid-Normal match?

- Identification of true somatic and true germline mutations
- Identification of fusions by DNA and RNA
- RNA-based gene expression on tissue
- High accuracy of TMB, MSI, and HRD with the patient's normal matched sample
- Accurate identification of reversal and Clonal hematopoiesis of indeterminate potential (CHIP) mutations
- High accuracy in detecting copy number variation (CNVs)
- Single comprehensive test to address complex molecular questions in cancer patients

OncoIndx[®] Prime⁺ in practice

- Challenges of managing treatment responses are addressed through tissue-liquid-normal match derived treatment choices.
- Tissue-liquid-normal match NGS report enhanced with AI and updated clinical studies help identify the most beneficial therapeutic options.

Clinical and Technical Validation of OncoIndx[®] Assay

Alteration	* PPV	* NPV	Accuracy	Specificity	Sensitivity
SNVs	100	100	100	100	100
Small INDELs	100	94.03	97.40	100	95.60
CNA	100	100	100	100	100
Fusions	100	96.43	98.48	100	97.44

Transcriptomics Validation Results


Sensitivity	* PPV	* NPV	Specificity	Concordance
86.7%	100%	100%	93.6%	95%

A total of 165 samples including reference standards, clinical, cross-laboratory and TCGA samples.

* NPV: Negative predictive value. PPV: Positive predictive value.

Designed for oncologists managing complex, advanced and refractory cancer patients.

- Proprietary panel of 1,000+ genes
- Tissue-liquid-normal matched
- RNA-based fusion, RNA expressed Variants, Splice Variants and gene expression (TBx)
- Longitudinal monitoring markers
- PD-L1 through IHC (TBx)
- CTC enumeration



When to use OncoIndx® Prime+ ?

01 Advanced Stage Cancers

Including progressive and treatment-resistant tumors

02 Rare Cancers

Sarcoma, Glioma to identify clinical-trial enrolment or off-label therapy

03 Cancers of Unknown Primary (CUP)

RNA-based gene expression analysis supports tissue-of-origin identification

04 Multiple or Dual Malignancies

Distinguish between independent primary tumors and metastatic outcome

05 Biomarker Profiling

Patient selection for treatment with Immunotherapy/ Targeted therapy


06 Ambiguous Histology

When pathology is inconclusive, multi-omic profiling refines molecular classification and treatment selection

Features of OncoIndx® Prime+ report

- In accordance with ACMG/AMP/ASCO/CAP guidelines
- Clear, unambiguous and actionable genomic findings.
- Based on highly reliable and internally curated databases.
- Precise reporting of tumor-derived, germline, reversal and CHIP mutations.
- ctDNA Tumor Fraction ensuring Lbx reliability

OncoIndx® Prime+ Workflow



OncoIndx® Prime+ is a well-validated Tissue-Liquid-Normal match NGS test to extract therapy information for the patient. It is possible on both tissue as well as blood sample.

1080 proprietary gene panel

This tissue-liquid- normal matched NGS test includes 1000+ genes of targetable/ therapeutic significance.

>99% on target coverage

The test has a mean coverage of 98-100% across all tested genes.

2000x / 10000x mean depth

The test is performed at a mean depth of 2000x for tissue samples and 10000x for liquid.



OncoIndx® TBx offers integrated genomic and transcriptomic profiling across solid tumors. Its unique design identifies key actionable insights for personalizing cancer treatment.

OncoIndx® TBx in practice

- Challenges of managing complex cancers are addressed through evidence-based treatment choice.
- NGS report enhanced with artificial intelligence and updated clinical studies help identify the most beneficial therapeutic options.
- An easy digital interface to interpret and understand complex cancer genomics.
- OncoIndx Tbx significantly improves the sensitivity and reliability of fusion detection in tumor samples

Where is OncoIndx® TBx needed?

- Cancers progressing on standard-of-care treatment.
- Cancers where personalised treatment could enhance survival.
- Cancers showing resistance to conventional treatment.
- Cancers where combination of target drugs and/or immunotherapy could benefit.
- Patients with rapidly progressing or aggressive tumors, where DNA only testing may miss actionable RNA fusions
- Cases where suspected gene fusions (ALK, ROS1, NTRK1/2/3 etc.) where RNA based detection increases sensitivity

Clinical and Technical Validation of OncoIndx® Assay

Alteration	*PPV	*NPV	Accuracy	Specificity	Sensitivity
SNVs	100	100	100	100	100
Small INDELs	100	94.03	97.40	100	95.60
CNA	100	100	100	100	100
Fusions	100	96.43	98.48	100	97.44

Validation summary for RNA sequencing Based Fusion Detection Assay


Sensitivity	PPV	NPV	Specificity	Precision
78.6%	100%	82%	100%	≥95%
Accuracy		LOD (empirical)		
≥95% Within validated input range		≥20ng RNA input with ≥60X coverage		

*PPV: Positive Predictive Value; NPV: Negative Predictive Value

Designed for oncologists managing complex, advanced and refractory cancer patients.

- Proprietary panel of 1,000+ genes
- Longitudinal monitoring markers
- PD-L1 through IHC (TBx)
- Traceability to ongoing clinical trials worldwide

OncoIndx® TBx



Key Actionable Genes Include:



Lung
EGFR, ALK, HRAS, NRAS, KRAS, BRAF, TERT, ROS1, MET, RET, NTRK1/2/3, ERBB2, PIK3CA



Breast
BRCA1, BRCA2, HER2, PIK3CA, MSH6, PMS2, NTRK1/2/3, MLH1, BRAF



Bladder
FGFR2/3, NTRK1/2/3



Colorectal
BRAF, KRAS, NRAS, HRAS, ERBB2, MLH1, MSH2, MSH6, PMS2, NTRK1/2/3



Gastrointestinal
ALK, HRAS, HER2, NRAS, KRAS, BRAF, TERT, ROS1, MET, RET, NTRK1/2/3



Melanoma
BRAF, KIT, ROS1, NTRK1/2/3



Prostate
BRCA1, BRCA2, MLH1, MSH2, MSH6, PMS2, AR, NTRK1/2/3



Thyroid
ALK, HRAS, NRAS, KRAS, BRAF, TERT, RET, NTRK1/2/3



Ovary
BRCA1, BRCA2, HER2, PMS2, MLH1, MSH2, MSH6, NTRK1/2/3

Features of OncoIndx® TBx report

- In accordance with ACMG/AMP/ASCO/CAP guidelines
- Clear, unambiguous and actionable genomic and transcriptomic findings.
- Based on highly reliable and internally curated databases.



Technical Specification

Specimen	Tissue
Panel Size	3 MB
Sample Type	1 FFPE tissue block with 20-30 % tumor content (not more than 3-4 months old)*
Sample Collection	FFPE tissue block
Sample volume/quantity	1 FFPE block with 20-30% tumor content
Input amount	DNA > 1ng RNA – 200pg – 1ug
Sequencing platform	Illumina
Read Length	2 x 150bp
Mean Sequencing Depth	2000x
Sequencing Chemistry	Hybrid capture
Indications	Solid Tumor
Analysis and Reporting	iCare™- AI based analysis platform with active reporting options

OncoIndx® TBx is an integrated NGS test that provides clarity to make impactful clinical decisions.

1080 proprietary gene panel

Comprehensive gene panel of targetable/therapeutic significance

>99% on target coverage

High quality, reliable results for therapy selection

2000x mean depth

Ultra deep sequencing for sensitive, reliable detection



OncoTarget® is cutting edge NGS Panel designed to identify targeted mutations in Bladder, Breast, Colon, GI, Lung, Melanoma, Ovary, Prostate, Rectal, and Thyroid.

OncoTarget® Test Outcomes

Genomic Alteration –


Easy to interpret - Tier classifications and appropriate therapeutic recommendations.

Note: PD-L1 expression on tissue is an add-on option if requested


Benefits of OncoTarget®

- Affordable genomic test with deep coverage of disease-specific guideline recommended genes.
- Easily comprehensible results.
- Improves accessibility of genomic tests.
- OncoTarget® is also **customizable** with ability to select **upto 20 genes** for each cancer type.

Types of alterations covered in OncoTarget®



How OncoTarget® tailors relevant information for patients?



Accelerate treatment decisions with OncoTarget®

- **Detect:** Clinically actionable guideline specific genomic mutations
- **Identify:** Relevant genomic targets and biomarkers (MSI)
- **Predict:** Disease prognosis and treatment responses

Key Actionable Genes Include:



Lung

EGFR, ALK, HRAS, NRAS, KRAS, BRAF, TERT, ROS1, MET, RET, NTRK1/2/3, ERBB2 (HER2), PIK3CA



Melanoma

BRAF, CD274, CDKN2A, NRAS, NTRK1, NTRK2, NTRK3, ROS1, MSI*



Prostate

BRCA1, BRCA2, CHEK2, MLH1, MSH2, MSH6, PALB2, PMS2, RAD51D, MSI*



Ovary

BRAF, BRCA1, BRCA2, MLH1, MSH2, MSH6, MTOR, NTRK1, NTRK2, NTRK3, RET, MSI*



Thyroid

ALK, BRAF, MLH1, MSH2, MSH6, NTRK1, NTRK2, NTRK3, PMS2, MSI*



Breast

BRCA1, BRCA2, ERBB2 (HER2), PIK3CA, MLH1, MSH6, PMS2, NTRK1/2/3, ESR1, BRAF, PALB2, RET



Bladder

MLH1, CDKN2A, ERBB2(HER2), FGFR1, FGFR2, FGFR3, NTRK1, NTRK2, NTRK3, PIK3CA, MSI, MSH6, NTRK1/2/3, PMS2, TP53, TSC1



Gastrointestinal

KIT, STK11, MLH1, MSH2, MSH6, PMS2, ALK, HRAS, ERBB2 (HER2), NRAS, KRAS, BRAF, ROS1, MET, RET, NTRK1/2/3, PDGFRA, NF1, FGFR fusions



Colorectal

BRAF, KRAS, NRAS, ERBB2 (HER2), MLH1, MSH2, MSH6, PMS2, NTRK1/2/3, TP53, EPCAM, PIK3CA, PTEN, RET, AKT1, BRAF, HRAS, MET, SMAD4

Technical Specifications

Specimen	OncoTarget® (LBx or TBx)
Sample Type	Plasma cDNA (BD Pax-Gene Tube) + 1 tube of 8.5 ml Peripheral Whole Blood [Note: ~8-10 days of washout period required if patient under therapy e.g., Chemotherapy, Radiation, Immunotherapy, etc.] OR 1 FFPE tissue block with 20-30 % tumor content. (Note: Not more than 3-4 months old)
Sample Collection	Peripheral blood using butterfly needles/ 22 gauge syringe
Sample volume/quantity	2 PAXGene tube + 1K2 EDTA peripheral whole blood OR 1 FFPE block with 20-30% tumor content
cfDNA/DNA input amount	> 1ng
Sequencing platform	Illumina
Read Length	2 x 150bp
Mean Sequencing Depth	2000x (tissue) / 10000x (liquid)
Sequencing Chemistry	Hybrid capture
Indications	Clinically diagnosed OncoTarget® associated cancers
Analysis and Reporting	iCare™- AI based analysis platform with active reporting options

Precision Targeted
Genomic Insight

Your Path to
Personalized Care



OncoHRD is a cutting-edge test that combines analysis of HRR gene alterations with genomic scar assessment to provide a comprehensive HRD status.



What is HRD (Homologous Recombination Deficiency)?

- HRD refers to a defect in the ability of cells to repair DNA double-strand breaks through homologous recombination.
- Commonly seen in cancers such as breast, ovarian, pancreatic, and prostate but is a pan-cancer aspect.

Why HRD: Some key concerns in ovarian cancer treatment?

- Homologous recombination deficiency (HRD) is identified in approximately 48% of ovarian cancer tumors, typically due to mutations specific to the tumor.
- Assessing HRD status is crucial, as it provides valuable insights into the potential benefits of PARP inhibitor therapy.





Figure. While 50% of High-grade serous ovarian carcinoma (HGSOC) samples have a positive HRD status, contributing genomic markers are found to be beyond BRCA mutations

Defects in the homologous recombination repair pathway

Ovarian Cancer

High-grade serous ovarian cancer (HGS-OvCa)

- 17% carried germline mutations in BRCA1/2
- 3% somatic mutations in BRCA1/2

Breast Cancer

Triple-negative breast cancer (TNBC)

- 20% had either somatic or germline mutations in BRCA1/2

Prostate Cancer

Metastatic castration resistant prostate cancer (mCRPC)

- ~19% had mutations in BRCA1/2, BRCA2, ATM, CDK12

Pancreatic Cancer

Pancreatic ductal adenocarcinoma (PDAC)

- 17% had germline mutations in BRCA2
- 24% germline or somatic mutations in BRCA1/2, PALB2
- 8% had mutations in ATM

Types of Genomic Scars included

Loss of heterozygosity (LOH)

Telomeric-allelic imbalance (TAI)

Large-scale state transitions (LST)

Advancing Precision Oncology with Next-Gen HRD Insights



Comprehensive Genomic Insights

Analyzes 50+ genes, detecting SNVs, INDELs, CNVs, and structural variants



Advanced HRD Scoring

Provides HRD scores with LOH, TAI, and LST for precise tumor profiling



AI-Powered Reporting

Utilizes iCARE™ for accurate, dynamic AI-powered reports




Optimized Treatment Guidance

Supports PARP inhibitors and platinum-based therapies for improved outcomes

Homologous recombination repair (HRR) gene

More than 50 HRR genes covered in OncoHRD® test

ARID1A	CDK12	FANCF	POLE	RPA1
ARID1B	CHEK1	FANCG	POLH	SMARCA1
ATM	CHEK2	FANCI	PPP2R1A	TP53
ATR	DNMT1	FANCL	PPP2R2A	TP53BP1
ATRX	ERCC1	FANCM	PTEN	WRN
AURKA	ERCC4	HDAC2	RAD50	XRCC1
BAP1	ERCC5	IRS1	RAD51	XRCC2
BARD1	EXO1	MLH1	RAD51B	XRCC3
BLM	FANCA	MRE11	RAD51C	
BRCA1	FANCC	NBN	RAD51D	
BRCA2	FANCD2	PALB2	RAD54B	
BRIP1	FANCE	POLD1	RAD54L	



HRD Biomarkers for
PARP Treatment Eligibility

HRD testing is applicable to a variety of cancers, including:



Breast Cancer

Especially useful for
BRCA-mutated or triple-
negative breast cancers.



Ovarian Cancer

A key predictor for
treatment with PARP
inhibitors.



Prostate Cancer

Helps identify those
who may benefit from
targeted treatments.



Pancreas Cancer

Guiding the use of
targeted therapies in case
of HRD positive cases.

Benefits of OncoHRD® test



Affordable Solution

Affordable NGS-based test with deep
coverage of 50+ genes, detecting
SNVs, INDELs, CNVs, and structural
variants



Swift Results

Easily comprehensible report generated with
the HRD score, genetic alterations using
AI-powered dynamic reporting through
iCARE™



Associated Drugs

Drug associations based on rigorous,
peer-reviewed literature ensuring
evidence-backed insights, linking treatments
to scientifically validated research.



Targeted Therapy

Targeted Therapy Selection by
identifying HRD-positive tumors,
this test enables oncologists to
select highly effective treatments
such as PARP inhibitors



Improved Survival

HRD-positive cancer patients often experience
better responses to targeted therapies, leading
to improved overall and progression-free
survival

When Is OncoHRD® Clinically Relevant for Your Patients?



Patients with a Family
History of Cancer




Patients with
Advanced Cancers



Patients Seeking
Personalized Therapy




OncoRisk® is a hereditary cancer test that represents a robust solution, examining a comprehensive panel of **74 genes** connected to a spectrum of **more than 11** distinct hereditary and familial cancers.



How can OncoRisk® predict the onset of cancer?

- Identify at-high-risk unaffected relatives of a cancer patient
- Prediction of risk which can lead to early detection
- Implement rigorous screening protocols for individuals at high risk and prognosticate the potential for cancer development.
- Timely intervention of cancer translates to enhanced treatment outcomes and prolonged periods of disease-free survival.
- Informed Surgical Choices, such as Prophylactic Mastectomy and Bilateral Salpingo-Oophorectomy



Designed for both oncologists and patients to uncover genetic mutations



More than 11 Distinct Hereditary and Familial Cancers



Pre and Post-test Genetic Counselling



AI-Powered Report with Actionable Insights



Key Actionable Genes Include:



Breast

ATM, BARD1, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, DICER1, MLH1, MRE11, MSH2, MSH6, MUTYH, NBN, PALB2, PMS2, POLD1, POLE, PTEN, RAD50, RAD51C, RAD51D, SMARCA4, STK11, TP53



Ovary

ATM, BARD1, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, EPCAM, MLH1, MSH2, MSH6, NBN, PALB2, PMS2, PTEN, RAD51C, RAD51D, TP53, STK11



Colorectal

APC, BMPR1A, CDH1, CHEK2, EPCAM, GREM1, MLH1, MSH2, MSH6, MUTYH, PMS2, POLD1, POLE, PTEN, SMAD4, STK11, TP53, AXIN2, NTHL1, MSH3



Gastrointestinal

APC, BMPR1A, CDH1, CTNNA1, EPCAM, MLH1, MSH2, MSH6, PMS2, SMAD4, STK11, TP53



Prostate

ATM, BRCA1, BRCA2, CHEK2, EPCAM, HOXB13, MLH1, MSH2, MSH6, NBN, PALB2, PMS2, PTEN, RAD51D, TP53



Pancreas

APC, ATM, BRCA1, BRCA2, BMPR1A, CDK4, CDKN2A, EPCAM, MLH1, MSH2, MSH6, PALB2, PMS2, SMAD4, STK11, TP53



Renal

BAP1, FH, FLCN, MET, PTEN, SDHA, SDHB, SDHC, SDHD, TSC1, TSC2, TP53, VHL



Skin

BAP1, BRCA2, CDK4, CDKN2A, EPCAM, FH, FLCN, MLH1, MIFT, MSH2, MSH6, PMS2, PTEN, TERT, TP53



Endocrine

MEN1, RET



Lung

EGFR, STK11



Benefits to Individuals Undergoing Germline Testing:


- Gain insight into your hereditary cancer risk, to take proactive steps for your health.
- Identify potential risks early, enabling timely interventions and personalised care.
- Make informed choices about your health, lifestyle, and medical decisions.
- Protect your loved ones by understanding the genetic basis of hereditary cancers.
- Receive a comprehensive report with actionable information for tailored care planning.
- Gain insights for informed risk-reducing surgeries that substantially lower hereditary cancer risks

Knowing Your Risk, Shaping Your Future

Your Genetic Shield Against Hereditary Cancers



OncoMonitor®MRD is a blood test designed for Longitudinal Cancer Monitoring and Minimal Residual Disease Detection, by measuring CTC & ctDNA.



OncoMonitor® TRM is a liquid biopsy-based solution designed for monitoring treatment response, emergence of resistance and tracking disease progression.

OncoMonitor®MRD Benefits



Minimal Residual Disease Detection

OncoMonitor's ability to identify trace amounts of ctDNA and CTCs allows for early detection of residual cancer cells post-treatment, enabling timely intervention and improved patient outcomes.



Prognostic Value

CTC and ctDNA levels have shown prognostic significance in various cancers. The test results can help predict disease progression and overall survival, aiding in informed decision-making for both patients and healthcare providers.

ctDNA epigenomics validation results

Sensitivity	90.5%
PPV	95%
Specificity	97.2%
NPV	94.6%
Concordance	94.7%
Limit of Detection (LOD)	0.05

Analysis done on 57 samples which includes Seraseq® reference standards and healthy individual samples

OncoMonitor® TRM Benefits



Longitudinal Monitoring

By regularly assessing CTC numbers and ctDNA load, oncologists can gain insights into the dynamics of the disease, its response to treatment, and potential disease relapse.




Treatment Efficacy Assessment

Changes in CTC counts and ctDNA levels can indicate whether the current treatment approach is successful or if modifications are necessary to enhance its impact.

For which cancer phase is OncoMonitor® most applicable?

Applicability Across Cancer Journey: Tracking Disease Progression and Relapse Detection Through Longitudinal Blood-Based Monitoring



Comparison of OncoMonitor® MRD and OncoMonitor® TRM

	OncoMonitor® MRD	OncoMonitor® TRM
Eligibility	Early-stage post-curative intent surgery	Advanced metastatic cancers
Features	Cellulomics and Epigenomics Captures CTC load & ctDNA load (Epigenomics)	Cellulomics, Epigenomics, and Genomics (mutation profiling) - Captures CTC load & ctDNA load (Epigenomics+Genomics) - Detects SNVs, INDEls, and CNVs - Covers 30 fusion genes with 40 corresponding partners on DNA alongwith MSI
Clinical Decisions	1. Escalation or de-escalation in adjuvant setting 2. Surveillance of tumor burden 3. MRD assessment	1. Early course correction 2. Predicting therapy response 3. Disease progression prediction
Test Type	Prognostic/Predictive test	Predictive test
Indications	Colon, Rectal, Breast, Bladder, Prostate, and Ovary cancers	All Solid Tumors

Test Performance Characteristics:

- **Gene coverage:** 99% mean coverage across all included genes
- **Sequencing depth:** 10000x mean depth
- **Variant Calling Databases:** Utilizes approximately 20 databases
- **DNA Methylation Profiling:** Covers over 3440 CpG sites
- **Sample Type & Volume:** Peripheral blood (2 PAXgene tubes, 1 K2EDTA tube; 10 ml each)
- **AI-Powered Reporting:** Enabled via iCARE™.



Key Actionable Genes Include in OncoMonitor TRM:



Lung

AKT1, ALK, BRAF, ERBB2 (Her2), EGFR, FGFR 1/2/3, HRAS, KRAS, MAP2K1, MET, NRAS, NTRK 1/2/3, PIK3CA, PTEN, RET, ROS1, TP53



Melanoma

BRAF, CTNNB1, GNA1, GNAQ, HRAS, KIT, KRAS, MAP2K1, NF1, NRAS, NTRK 1/2/3, PDGFRA, PIK3CA, PTEN, TP53



Colorectal

AKT1, BRAF, HRAS, KRAS, MET, MLH1, MSH2, MSH6, NRAS, NTRK 1/2/3, PIK3CA, PMS2, PTEN, SMAD4, TP53



Gastrointestinal

BRAF, KIT, KRAS, MET, MLH1, MSH2, MSH6, NTRK 1/2/3, PDGFRA, PMS2, SMAD4, TP53



Breast

AKT1, AR, BRCA1, BRCA2, ERBB2 (Her2), ESR1, FGFR 1/2/3, NTRK 1/2/3, PALB2, PIK3CA, PTEN, TP53



Bladder

MLH1, MSH2, MSH6, NTRK 1/2/3, PMS2, TP53, TSC1



Ovary

BRAF, BRCA1, BRCA2, HRAS, KRAS, NRAS, NTRK 1/2/3, PALB2, PDGFRA, TP53

**Pioneering
Precision Monitoring
Navigating Cancer's
Ebb and Flow**

Longitudinal Monitoring

Assess disease
dynamics
over time

Minimal Residual Disease detection

Uncover residual
cancer cells
post-treatment

AI powered Real-time results

Swift and accurate AI
powered comprehensive
reporting



1Cell.Ai™
Ai-Powered Precision Oncology



Learn More at 1Cell.Ai
Mail - info@1cell.ai | Tel - 1800 123 12355
1Cell.Ai © Copyright 2025. All rights reserved | Scan QR to know more





OncoIndx[®] 360[°] Endometrium is the first upfront, integrated molecular test that delivers a single, biology-driven snapshot with clear, guideline aligned treatment direction in endometrial cancer.


One test. One tissue sample. One integrated molecular snapshot:

International guidelines (ESGO/ESTRO/ESP) and clinical trials now clearly demonstrate that molecular classification outperforms histology alone in predicting outcomes and guiding therapy in EC. Yet, today, critical molecular tests are ordered sequentially across multiple labs—delaying decisions, exhausting tissue, and leaving key categories like NSMP poorly interpreted.

Before:

Current Practice


Multiple tests | Multiple labs | Multiple delays



After:

OncoIndx360[®]

Integrated | Upfront | Actionable



From Ambiguity to Standardized Decisions:

- Receive complete molecular risk stratification upfront, before adjuvant planning
- Avoid overtreatment in POLE-mutated disease despite high-grade histology
- Reduce repeat biopsies and tissue exhaustion
- Confidently tailor adjuvant radiation decisions based on true recurrence risk
- Integrate molecular risk into escalation or de-escalation strategies

Empowering oncologists to decode complex, advanced, and refractory endometrial cancers through integrated molecular insight

- Upfront Integrated Profiling (NGS + IHC)
- Single Unified Report
- Biology-Driven Interpretation
- WHO aligned treatment direction

Selected Actionable Molecular Classifications Include:

Molecular Subtype	Prognosis	Clinical Implication
POLE-Ultramutated	Excellent	Safe treatment de-escalation; avoid unnecessary chemotherapy despite high-grade histology
MMRd / MSI-High	Intermediate (therapy-responsive)	Prioritize immunotherapy; moderate recurrence risk without IO
TP53-Abnormal (p53abn)	Poor	Treatment escalation; intensive chemotherapy ± HER2-targeted therapy
NSMP – ER Positive	Favorable-Intermediate	Hormone-based therapy; reduced reliance on cytotoxic chemotherapy
NSMP – ER Negative	Intermediate–Poor	Standard cytotoxic therapy; closer surveillance
HER2-Positive	Poor (improves with targeted therapy)	Precision escalation with HER2-targeted agents


Technical Specifications:



Sample type: One FFPE Tissue block



Indication : Endometrial cancer




Technology: NGS based integrating profile



Turn around time : 12-14 days



LOD (Limit of Detection) – 1% VAF



Output: Molecular classification, risk stratification, therapy guidance



Analysis and Reporting : iCare™



Clinically validated : With high sensitivity

Features of OncoIdx® 360° Endometrial Report:

- Integrated NGS, IHC results in a single, unified report, in addition to Biologically Relevant Biomarker Findings (RUO), offering research insights.
- Comprehensive actionable WHO-aligned biomarker coverage (POLE, TP53, MMR, ER, HER2)
- Clear prognosis and therapy guidance including immunotherapy, targeted, hormone therapy options.

OncoIdx® 360° Endometrium is an integrated NGS and IHC test, that provides clarity to make impactful clinical decisions

Proprietary gene panel

Comprehensive gene panel of targetable/therapeutic significance

>99% on-target coverage

High quality, reliable results for therapy selection

2000x mean depth

Ultra deep sequencing for sensitive, reliable detection

